

## ColoQuick Pasteur + MAXX4

# MENU CODES

<b>PROGRAM SETUP</b>	
1110, 1111, 1112	Pasteurization setup
1120, 1121, 1122	Thaw setup
1130	Keep warm setup
1140, 1141, 1152	Cool setup
2222	User setup
2223	Hot water setup
2240	Calibration
2242	Calibration status
2999	Clear Error Queue
3000	Show log
3001	Show Error Queue
3010	Clear log
7500	Load USB settings
7501	Save USB settings

<b>FACTORY SETUP</b>	
8191	Factory test
8192	Set product ID
8195	Enable motor PWM

<b>MACHINE SETUP/TEST</b>	
1280	Wear status
9990	Test mode
9991	Sensor configuration
9992	Enable/disable func
9999	Software upgrade

